



Commercial Products    **Military Products**    Tips & Uses    Reviews    Reports    Where to Buy    What's New    Contact Us

## Military Products

### XiniX Military KlearWater Multi-Use Water Treatment

XiniX's Military version of KlearWater is formulated for multi-use and is the perfect companion for the field, for water, hands, surfaces—all cleaning needs. Available in bulk pack with or without snap container.

#### Instructions on label:

Add ½ cap per liter (quart) of water and wait 15 minutes  
Add 1 cap per liter (quart) of water if Cloudy or suspected fouled by animal waste and wait 30 minutes  
Use to wash/rinse hands, fruits and vegetables, utensils  
Field expedient for sores, cuts, scrapes, blisters, herpes, boils

**XiniXtremedic**—Antiseptic Formula  
50 ml glass bottle with cap

#### Instructions on label:

Dilute solution 3 to 1 with clean water to Flush and irrigate wounds.  
For open Wounds, blisters, cuts and scratches: Clean with XiniXtremedic and cover with a sterile dressing.  
Use for cleaning and cold sterilization of surgical instruments in the absence of other cleaning methods.  
Spray skin prior to suturing.

### XiniX SprayDose

Sprayer packed in plastic tube with cap  
Use with 50 ml Xtremedic

- ACTIVE CHLORINE DIOXIDE IN DEIONIZED WATER

#### Instructions on label:

Squirt a Nalgene bottle with a 3 to 1 dilute solution to wash wounds, sores, scratches  
Spray full strength prior to suturing  
Spray full strength after flushing, cover with dressing  
Field expedient for cold sterilization

### XiniXtreme

Water disinfectant

- Disinfects Water FAST
- Maintains Bladder System while disinfecting
- Handy single-use bottle
- Field expedient for hand wash and other field needs
- Perfect Companion for 3 liter Bladder
- Disinfects and Maintains
- Error-Proof Dosage, Single-Use
- Fastest/Lightest/Most Effective.
- Field Tested 16 months by ASOC
- Exclusive recommendation of JSOC Medics

#### Instructions on label:

Instill bottle (5 ml) contents into 3 liters  
Water may be turbid but not dark  
Mix and allow 15 minutes to neutralize pathogens  
Double time and dose if animals present  
Field expedient for wounds and washing  
Active when yellow, keep from eyes  
Do not drink concentrate

Same formula used by professional extreme athletes worldwide! Available in a box of 104.

### Xinix Surface Concentrate with Reusable Double Acting Sprayer/Mister

Cleans all surfaces more powerfully than bleach

- Deploy with Xinix Surface Concentrate, to produce XiniXtreme cleaner in the field with potable water
- (at least) 5X more powerful than bleach for cleaning in the field
- Compact for field kitchen, field medic needs, latrine

Available in bulk pack or 24 packs with separators for rough handling.

**Xinix KlearWater products** are fast acting, lightweight, multi-use, economical, and don't add unpleasant taste or smell to water. Please contact our [sales department](#) for inquiries or to establish an account. We do not repackage and



